ASPEL S.A.

32-080 ZABIERZÓW, os. Sienkiewicza 33

tel. (012) 285-22-22, fax (012) 285-30-30

www.aspel.com.pl

Elektrokardiograf

AsCARD A4

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WYDANIE IX

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tel. +48 12 / 285 22 22, fax +48 12 / 285 30 30

www.aspel.com.pl

**Operation Manual
Electrocardiograph
AsCARD A4**

Zabierzów, POLAND

October 2004

EDITION IX

CONGRATULATIONS

- We are pleased that you have chosen the AsCARD A4 electrocardiograph;
- Choosing the services offered by ASPEL you have proved that you do not approve compromise and you want the highest quality.
- ASPEL has the pleasure to present the AsCARD A4 electrocardiograph which is a result of years of surveys and experiences gathered through close contact with our customers. You have chosen high durability and efficiency as these are the features of AsCARD A4.
- Read this Manual carefully as it includes guidelines concerning the safe installation, use and maintenance and several practical instructions how to optimise this apparatus application.
- Store this Manual for reference.

1. APPLICATION

The electrocardiograph AsCARD A4 is a high quality unit, three, six and twelve channels ECG unit which enables to make electrocardiogram in full 12 leads. The ECG unit maintenance is very simple and the modern keyboard makes the operation much easier. The record is made in manual or automatic mode with the possibility to make analysis and interpretation. The unit can be supplied from the net with 230V voltage or from internal accumulators built in the unit. The accumulator capacity makes possible to print about 60 examination in auto mode. Low magnitude and small size are the additional features which enable easy transport. It is possible to carry out spirometry examinations and to send ECG examination by fax with the use of this ECG unit. The electrocardiograph is designed for carrying out ECG examinations at the patients home and in all health-care institutions.

1. The person operating the electrocardiograph ought to read the Operation Manual carefully before starting the unit operation.
2. The Operation Manual will help the user to operate and maintain electrocardiograph correctly.
3. The efficient functioning and fulfilling the tasks the user reaches after obeying the Operation Manual descriptions.
4. It is forbidden to use the unit with a damaged patient cable.
5. Every repair ought to be made by Authorised Services.
6. The manufacturer is not responsible for any damages resulting from not obeying the Operation Manual.
7. If the electrocardiograph is used in combination with other equipment assure that:
 - all connections are in compliance with international IEC safety standards.
 - all connected units are in compliance with IEC standardsIf the above standards are not obey it can cause the safety hazard for which the manufacturer is not responsible.
8. Electrocardiograph used together with a cardiostimulator or other electric stimulator do not cause any safety hazard of the patient or the user.

9. The efficiency of the unit and accessories shall be periodically controlled (min. once a year). If there are any irregularities in functioning of the unit, one shall contact Authorised Service, every time when they appear.
10. For the purpose of long-term archivisation it is requested to make copy of printed examination (photocopy, scan). Thermoactive paper is too sensitive for the variable environmental conditions. It can have a result in illegibility of printout kept for a long period.


1. 1 Description of unit.

The device is made up of modern microprocessor technology and all parts are produced by reliable international companies. The unit has a thermal printer, equipped with a head of high resolution to print the ECG leads.

The cover made out of ABS resins and the membrane keyboard facilitates cleaning.

1.2 Picture of unit.

AsCARD A4

- 1. - paper box
- 2. - RS232 connection
- 3. - power supply socket
- 4. - switcher on/off
- 5. - LCD display
- 6. - keyboard
- 7. - patient cable socket

1.3 Basic Technical Parameters

Parameter	AsCARD A4			
	Model A	Model B	Model C	Model D
Dimensions: (LxWxH) [mm]	310x230x66	310x230x66	310x230x66	310x230x66
Weight : ca [kg]	2,6	2,6	2,7	2,7
Power supply	230 V $\pm$ 10% 50 Hz $\pm$ 2%	230 V $\pm$ 10% 50 Hz $\pm$ 2%	230 V $\pm$ 10% 50 Hz $\pm$ 2%	230 V $\pm$ 10% 50 Hz $\pm$ 2%
Power consumption: [VA]	max 20	max 20	max 20	max 20
Internal accumulator *HITACHI HP1,2-12	12V 1,2 Ah	12V 1,2 Ah	12V 1,2 Ah	12V 1,2 Ah
ECG leads	12	12	12	12
Sensitivity : [mm/mV]	2,5/5/10/20	2,5/5/10/20	2,5/5/10/20	2,5/5/10/20
Record speed : [mm/s]	5/10/25/50	5/10/25/50	5/10/25/50	5/10/25/50
Paper width : [mm]	112	112	112	112
Vertical record definition : [point/mm]	8	8	8	8
Impedance [M Ω]	>10	>10	>10	>10
CMRR [dB]	>100	>100	>100	>100
Frequency range: [Hz]	0.05 $\pm$ 150	0.05 $\pm$ 150	0.05 $\pm$ 150	0.05 $\pm$ 150
Digital filters : [Hz]	50 and 25, 35	50 and 25, 35	50 and 25, 35	50 and 25, 35
LCD Display	2x24 marks	2x24 marks	graphical 240x128	graphical 240x128
Keyboard	function	alphanumeric	function	alphanumeric

*the accumulator is changeable by Authorised Service only.

Using safety:

- security class : I
- security degree :CF(EN 60601-1)

ECG input is secured against the defibrillating impulse.

Fuses:

- | | | | |
|---|---------------|---------------|----------------------|
| • | mains line | T500mA/L/250V | <i>changeable by</i> |
| • | mains neutral | T500mA/L/250V | <i>the user</i> |

at the back side of the unit cover

- | | | | |
|---|---------------------|--------------|--------------------------------|
| • | battery charging B1 | T6,3A/L/250V | <i>changeable by</i> |
| | | | <i>Authorised Service only</i> |

Attention !

The electrocardiograph AsCARD A4 is the unit belonging to the A class. If used in households the unit may cause radio-waves disturbance. If the disturbance occurs the user may be forced to take precautions and bear responsibility for them.

Electrocardiograph should be connected to single-phase electrical system with ground connection.

1. 4. Manufacturer

ASPEL S. A.
os. H. Sienkiewicza 33
32-080 Zabierzów / POLAND
tel. + 48 12 / 285 22 22, fax. + 48 12 / 285 30 3

1. 5 Versions of unit.

AsCARD A4 is produced in four models A, B, C and D:

Software version:

- * The measurements of complexes and waves' time
- * Interpretation (includes measurements)

OPTIONS:

I – interpretation

IR – interpretation and RS connection

The software versions are marked in the unit number.

The marking consists of:

--	--	--	--

Factory number

--	--

Year of production

--	--

Option

2. WORKING CONDITIONS

Working conditions of electrocardiograph AsCARD A4

- temperature +5 to + 40 °C

- relative humidity < 95%

- the air should not be polluted with corrosion ingredients.

3. FEATURES

The electrocardiograph AsCARD A4 allows the recording of 12 standard ECG leads.

The following modes are possible:

- **3-channel manual record**

this mode allows the recording of selected ECG lead. It is possible to set required record speed, sensitivity, filters on/off: 50 Hz (power supply disturbances), 25 Hz, 35 Hz (muscular disturbances). It is possible to print one of four ECG leads I-II-III, aVR-aVL-aVF, V1-V2-V3, V4-V5-V6.

- **6-channel manual record**

this mode allows the recording of selected group of leads. It is possible to set required record speed, sensitivity, filters on/off: 50 Hz (power supply disturbances), 25 Hz, 35 Hz (muscular disturbances). It is possible to record one of two groups of leads - limb : I-II-III-aVR-aVL-aVF and precordial: V1-V2-V3-V4-V5-V6.

- **12-channel manual record**

this mode allows to record 12 leads I-II-III-aVR-aVL-aVF-V1-V2-V3-V4-V5-V6. It is possible to set required record speed, sensitivity, filters: 50 Hz (power supply disturbances), 25 Hz, 35 Hz (muscular disturbances).

- **automatic record**

this mode allows the recording of 10 sec of ECG signal from all 12 leads in three, six channel. It then allows an automatic analysis including time intervals, segment amplitude measurement and the calculation of the electric axes. The full report printout includes ECG signal record, averaged P-QRS-T complexes, measurements and calculations results. Record parameters are the same as in manual record.

- **spirometry examination record**

this mode allows the recording of the spirometry examination including Flow Volume Capacity.

- **copy of automatic record**

this mode allows the recording of copy of automatic examination stored in the unit memory.


It is possible to print a few dozen examples of the examination on 112 mm.

3. 1 How to prepare AsCARD A4 to work

The mains operation:

Connect the electrocardiograph AsCARD A4 to the power source. Put the main cable to the mains socket (3) placed on the back side the unit and to the power source 230V. The unit is switched on/off with the use of the switcher (1) placed on the back side of the unit.

To switch the unit on press the  key. The display shows:

AsCARD ver 3.xx


Internal accumulator operation.

If the unit hasn't been connected to the mains supply it works on internal accumulator. The unit switches off after 3 min, if any key has't been pressed. To switch the unit on

again press the  key.

Note !

The internal accumulator starts to charge automatically once the unit has been connected to the mains and switched on with the use of on/off switcher.


View of the back side of the unit

1. The mains switcher
2. Fuse
3. The mains socket
4. RS socket

3. 2 How to put ECG paper into the unit.

To put the paper into the unit, remove the cover of the paper box then press the unit power supply. "OUT OF PAPER" appears on the display. Put the roll inside the box and slip the paper through the driving roller. The unit should take the paper automatically. Then put the cover of the paper box on.

To prevent the linear head from pollution and the unit from blocking, it is necessary to use thermoactive and dustfree paper, which is always available at Aspel.

Technical Specification of ECG paper:

- Kind of paper	thermoactive, free of wax
- Kind and colour of overprint	millimetre net, orange
- Width of paper	112 mm $\pm$ 1 mm
- Diameter of a roll of paper	45 mm $\pm$ 1 mm
- External diameter of the paper shaft	16 mm $\pm$ 0,5 mm

Note !

The producer of electrocardiographs do not take the responsibility for any defects caused by using the other kind of paper.

Power supply signalisation :

a) the mains

if the unit is connected to the net supply it is signalised on the display with "→"

RECORD	I	II	III
25mm/s	10mm/mV	→	

b) accumulator

if the unit works on an internal accumulator, the switch on and start of auto and manual examination are signalised with acoustic "beep". If the "⎓" sign is displayed it signalises the accumulator discharge. The accumulator has to be charged for 12 hours.

ATTENTION !


There is a battery containing the lead in the apparatus. Don't throw the used up battery to a rubbish - bin! Subject it to recycle according to the obligatory regulations.

3. 3 Arrangement of the keys.

The arrangement of the keys is presented below:


Model A and C

Function of the keys

-  power supply ON key
-  power supply OFF key
-  automatic examination start key
-  manual examination start key
-  stop paper key
-  rhythm key
-  choosing lead key
-  record sensitivity key
-  copy key
-  filters key
-  record speed key
-  registration mode choosing key
-  spirometry key
-  fax transmission key

Model B and D

Model B and D have the additional keys

filter key of muscular disturbances

filter key of the mains disturbances

1mV pattern printout key


Menu - the configuration mode key

There are two kinds of **LCD display**:


Model A and B:

(2 lines x 24 marks) presents all the information needed in order to operate the unit.
The example display is shown below:

Model C and D:

Graphical LCD display presents 3 leads (I, II, III) and two text line with information needed for unit operation. The example display is shown below:

3. 4 Arrangements of the electrodes.


The electrocardiograph AsCARD A4 is equipped with a 10-electrode patient's cable. In order to receive 12 standard leads of Einthoven, Goldberger, Wilson, the electrodes should be attached in the following way:

Limb electrodes:

R - red	- right arm
L - yellow	- left arm
F - green	- left leg
N - black	- right leg

Precardiac electrodes :

C1; white-red	the fourth intercostal space on the right side of the sternum
C2; white-yellow	the fourth intercostal space on the left edge of the sternum
C3; white-green	half way between C2 and C4
C4; white-brown	the fifth intercostal space in the left central clavicular line
C5; white-black	straight line from the C4 point perpendicularly to the left front auxiliary line in the intersection point with this line.
C6; white-violet	at the same level as C5, but in left central auxiliary line.

How to clean the electrodes.

To clean the electrodes use a duster moistened with water and soap.

If the patient has been simultaneously connected to the several units, one should

remember to estimate any possible risks that may have the result from leakage current summing for each of used units.

Pay attention !

During the electrodes connection, the conducting parts of electrodes and patient cable should not be in contact with each other or with other metallic parts, including ground connection. The electrocardiogram shall appear in no longer than 10 seconds after de-fibrillator impulse.

4. EXAMINATION PERFORMANCE

4. 1 How to perform manual examination.

Switch on the unit with the use of the  key.

The display shows:

STOP	I	II	III
25mm/s	10mm/mV		

The unit supervises the electrodes contact with the patient's body and signalises poor contact with **INOP**.

INOP	I	II	III
25mm/s	10mm/mV		

When **INOP** is signalised, the ECG leads are printed as straight lines. It doesn't matter which electrode is badly connected. Check the contact on each electrode. When done:

Choose the record speed

Choose the record sensitivity

Choose the lead to record

Choose the registration mode

Press the  key. The unit starts electrocardiogram registration - the chosen ECG leads, speed, sensitivity and filters are presented on the display.

PRINT	I	II	III
25mm/s	10mm/mV		

During the electrocardiogram registration it is possible to change such parameters as:

- ECG leads
- record speed
- record sensitivity
- filters switching on or off
- record mode

Every parameter's change is automatically described on the display and registered on the electrocardiogram.

ECG registration lasts until the  key is pressed - then registration is broken off.

Press the  key to switch off the unit.

NOTE!

The unit switches off automatically. If the unit is battery operated and no keys are pressed or nothing is printed, it switches off after 3 minutes.

To switch the unit on press the  key.

4. 2 How to perform automatic examination.

Switch on the unit with the use of the  key. The display shows:

STOP	I	II	III
25mm/s	10mm/mV		

The unit supervises the electrodes contact with patients body and signals poor contact with **INOP**,

INOP	I	II	III
25mm/s	10mm/mV		

When **INOP** is signalised, the ECG leads are printed as straight lines. It does not matter which electrode is badly contacted. If the INOP is displayed the auto examination will not be carried out. Check the electrodes contact before the examination start.

Choose the record speed (25, 50 mm/s can only be chosen)

Choose the record sensitivity

Choose the lead to record

Choose the registration mode

Press the  key. The display shows:

AUTO	I	II	III
25mm/s	10mm/mV		

The unit stores the ECG in the memory simultaneously from all 12 leads. After memorising 10 sec course, the unit starts the ECG printout (described before the examination start). The length of ECG leads is described in the unit configuration.

The examination can be interrupted by pressing the  key.

If in configuration parameters the auto examination with analysis was set or during data collecting the  key was pressed, the unit starts to calculate the ECG parameters. The display shows:

PLEASE WAIT
*** ECG ANALYSIS ***

After 20 sec. the unit starts to print the report including:

- P-QRS-T complex of the most numerous class with marking beginnings and ends of segments. The lower line shows the class size. The record speed is always 50 mm/s.
- table with the values of segments' amplitude in all leads
- written interpretation of P-QRS-T complex's calculated parameters.

After finishing the printout the unit goes back to the beginning position:

STOP	I	II	III
25mm/s	10mm/mV		

NOTE!

The measurement results are reliable if on the QRS class representative record the beginnings and ends of segments are marked correctly and the number of detected complexes is >2.

Press the  key to switch off the unit.

4. 3 How to perform copy of automatic examination

Switch on the unit with the use of  key. The display shows:

STOP	I	II	III
25mm/s	10mm/mV		

Choose the record speed

Choose the record sensitivity

Choose the registration mode

Press the  key. If there is any correct auto examination in the memory the display shows:

EXAMINATION	NR 5
10:20:00	00.12.15

The unit memorises 10 last automatic examinations.

If there is the examination in the unit memory, the time and date are displayed in the second line. To change the number of the examination use of "↓, ↑".keys.

To print memorised examination simply press the  key. The display shows:

WAIT	NR 5
10:20:00	00.12.15

The unit starts to print the memorised electrocardiogram in the mode described before the print start. The duration of ECG leads is described in the unit configuration. Examination can be interrupted by pressing the  key.

If in configuration parameters the auto examination with analysis was set, then the unit starts to print its results. When the printout is finished the unit comes back to the start position. The display shows:

EXAMINATION	NR 5
10:20:00	00.12.15

To escape from copy mode, press the  key. The display show:

STOP	I	II	III
25mm/s	10mm/mV		

Press the  key to switch off the unit.

4. 4 How to record the rhythm

Switch on the unit with the use of  key. The display shows:

STOP	I	II	III
25mm/s	10mm/mV		

Press the  key. The unit goes into the rhythm mode. The display shows:

STOP	HR=60		
25 mm/s	10 mm/mV		

Choose the record speed

Choose the record sensitivity

Press  key. The unit starts the registration of I, II ECG leads and the values of heart rate frequency HR.

The registration lasts till the  key has been pressed.

To escape from rhythm mode, press the  key. The display show:

STOP	I	II	III
25mm/s	10mm/mV		

Press the  key to switch off the unit.

5. HOW TO CHANGE REGISTRATION PARAMETERS

Switch on the unit with the use of the  key.

The display shows the parameters set before switching off the unit.

STOP	I	II	III
25mm/s	10mm/mV		

Information showing on the display is printed on the paper as description of registered

ECG leads. Every parameters change, during the registration, causes the print of ECG status on the printout.

The change of sensitivity, leads or filters causes additional printout of 1 mV pattern.

5.1 Record speed.

It is possible to record in 5; 10; 25; 50 mm/s. The chosen value is presented on the display. To change the speed use the  key.

5.2 Record sensitivity.

It is possible to set 2,5; 5; 10; 20 mm/mV record sensitivity. The chosen value is presented on the display. To change the sensitivity use the  key. Every change of the sensitivity, during the registration, causes the printout of 1mV pattern and present status.

5.3 ECG leads.

To change the leads use the  key. The leads are changing from I to V6. Every change of leads, during the registration, causes the 1mV pattern printout with the leads description.

5.4 Additional filter.

It is possible to set additional ECG signal:

50 Hz	mains disturbances filter
35 Hz	muscular disturbances filter
35 Hz/50 Hz	muscular and the mains disturbances filter
25 Hz	muscular disturbances filter
25 Hz/50 Hz	muscular and the mains disturbances filter

It is possible to choose 60 Hz instead of 50 Hz. It can be made in "Configuration mode" - "Mains filter" (Chapter 6.).

Model A and C:

To change filter parameters use the  key.

Model B and D:

To switch 25, 35 Hz filter press the  key, to switch 50 Hz press the  key.

5. 5 Registration mode.

It is possible to make ECG examination in below given modes:

3 channels sequenced printout of 3 leads

6 channels printout of 6 leads

12 channels printout of 12 leads

To change the modes use the  key.

3-channel

STOP	I	II	III
25mm/s	10mm/mV		

6-channel

STOP	I	II	III	aVR	aVL	aVF
25mm/s	10mm/mV					

12 channel

STOP	12 LEADS
25mm/s	10mm/mV

6. HOW TO CONFIGURE THE UNIT PARAMETERS.

The electrocardiograph AsCARD A4 has a battery undercutting memory where is written: the data describing the starting position of the unit (after switching the unit on), the last automatic examination and report format.

The parameters relating to the configuration printout, isoline filtration and ECG examination type can be changed in configuration mode only.

The meaning of the keys are changed in configuration mode. Arrows placed on the left side of the keys are used.

- to choose the option "←, →"
- to confirm "↓"
- to cancel and exit from configuration mode "↑".

Switch off the unit with the  key.

To enter the configuration mode:

Model A and C:

Press the  key and the  key at the same time.

Model B and D:

Press the  key The display shows:

CHOOSE LANGUAGE:
>OK ENGLISH

It is possible to set one of the languages given below.

The list of languages:

- | | |
|------------------|-----------------|
| 1. ENGLISH | 2. POLSKI |
| 3. DEUSTCHE | 4. FRANCAIS |
| 5. ITALIANO | 6. ESPANOL |
| 7. HRVATSKI | 8. TURKCE |
| 9. CESKY | 10. ROMANA |
| 11. SLOVENSKY | 12. LATVIESU |

To set the language use "→" key. Choose the required one using "↓, ↑".keys, set ">" mark using "←" key in OK position and accept it with the use of "↓" key.

Next parameter appears in **Model B** and **D** (alphanumeric keyboards):

ENTER SURGERY NAME

You can insert the "Surgery name" with the use of alphanumeric keyboard. The description is printed once automatic examination is performed.

After pressing "↓" key, the display shows:

INOP CONTROL
>FULL SIMPLIFIED

It is possible to chose the mode of electrodes control connection: FULL or SIMPLIFIED (software ver. 1.43 and next).

The "FULL INOP CONTROL" - the unit automatically controls the isoelectric line position. If any of electrode is badly connected, the ECG recording is impossible. After correcting the connection, the electrocardiogram appears within 3s.

The SIMPLIFIED option is used in the case if you want to carry out the ECG examination without connection of all the electrodes (e.g. only limb electrodes are connected).

The "SIMPLIFIED INOP CONTROL" is possible only in manual record. In this option the electrocardiograph controls the connection of the electrodes and signalises the lack of connection with INOP. But the filters which controls the isoelectric line position are not switched on. In this setting electrocardiogram appears within 15s since the last electrode has been connected. In model B it is possible to shorten that time

with the use of the  key.

To set the option use " $\leftarrow$, $\rightarrow$ " keys. Set ">" mark in required position and accept with " $\downarrow$ " key, then next parameter is shown.

AUTO MODE WITH ANALYSIS:

>YES NO

If **YES** is chosen, then after all ECG leads that have been printed out in the automatic examination, the unit will start the analysis of the memorised leads and the calculated results will be printed.

If **NO** is chosen in automatic examination mode, the ECG will be printed and the unit will return to the start position.

To set the option use " $\leftarrow$, $\rightarrow$ " keys. Set ">" mark in required position and accept with " $\downarrow$ " key.

Next parameter appears on the display:

ECG INTERPRETATION:

>YES NO

Option **YES** - after calculation and table printout the interpretation is printed. To accept the choice press " $\downarrow$ ", the unit goes to the next parameters.

To set the option use " $\leftarrow$, $\rightarrow$ " keys. Set ">" mark in required position and accept with " $\downarrow$ " key.

Model A and C:

AUTO MODE PRINTOUT

>OK 5 [s]

To set the option use " $\rightarrow$ " key. Choose the time using " $\downarrow$ $\uparrow$ " keys. Set ">" mark using " $\leftarrow$ " key in OK position and accept with " $\downarrow$ " key.

Model B and D:

AUTO MODE PRINTOUT
5 [s]

Insert the time from the alphanumeric keyboard and press "↓" key. The display shows:

ISOLINE FILTRATION:
> YES NO

Option **YES** causes ECG leads filtration.

Option **NO** switches off the filtration.

To set the option use "←, →" keys. Set ">" mark in required position and accept with "↓" key, then next parameter is shown.

MAINS FILTER
> 50 Hz 60 Hz

To set the option use "←, →" keys. Set ">" mark in required position and accept with "↓" key, then next parameter is shown.

PRINTOUT PRECISSION
> 1 2 3

To set the option use "←, →" keys. Set ">" mark in required position and accept with "↓" key, then next parameter is shown.

Model A and C:

TIME HOUR MIN
OK> 09 45

To set the time use "←, →" keys. Place ">" mark in hour then minute position. Set the required time with the use of "↓ ↑" keys. Then set ">" mark in OK position and accept with "↓" key.

Model B and D:

TIME 09 : 45

To set the time use the keys on the alphanumeric keyboard. Set hour, press "space" then set minutes. To accept the choice press "↓".
The display shows:

Model A and C:

DATE	YEAR	MONTH	DAY
>OK	98	12	11

Model B and D:

DATE : 98 : 12 : 11

The date is set in the same way as time. Press "↓" key. The display shows next parameter:

FAX MODEM CLASS
> 2 2.0 NOKIA OFF

Setting the fax modem class allows to adjust the format of sending data to the mode that fax modem uses. The „OFF" position switches the modem option off.
To set the option use "←, →" keys. Set ">" mark in required position and accept with "↓" key. The display shows:

STOP	I	II	III
25mm/s	10mm/mV		

7. HOW TO TRANSMIT ECG EXAMINATION BY FAX

Switch on the unit with the use of  key.

STOP	I	II	III
25mm/s	10mm/mV		

To send ECG examination by fax you have to prepare:

- AsCARD unit with RS modem cable (produced by Aspel),
- modem that works like fax in class 2 or 2.0,
- telephone connected to the telephone line.

For the purpose of sending ECG data, the fax modem has to have the possibility to work

in class 2 or 2.0. Check the class of the fax modem you use in its operation manual. Set the same class in AsCARD unit - see "Configuration mode", chapter 6. Connect the fax modem to the electrocardiograph with the use of RS modem cable and to the telephone line. Make sure that there is an auto examination saved in the unit.

Press the  key. The display shows:

MODEM
INITIALISATION

Model A and C

After modem initialisation the display shows:

AFTER FAX SIGNAL
PRESS COPY

Choose the fax number with the use of telephone keyboard. If there is fax signal given - press the  key.

Model B and D

(for the "receiving fax" working in auto mode)

After modem initialisation the display shows:

FAX :

Choose the fax number with the use of alphanumeric keyboard of ECG unit. Accept with Enter "↓". If there is fax signal given - the copy of examination is sending automatically.

Model B and D

(for the "receiving fax" working in manual mode)

After modem initialisation the display shows:

FAX :

Press Enter "↓" and act like in Model A and C.

Next operations are the same for both models A and B.

CONNECTING

Sending ECG examination:

ECG TRANSMISSION

If the transmission is finished, the display shows:

TRANSMISSION OK

The transmission can be interrupted at any time by pressing the  key. If there is an error the display shows:

TRANSMISSION ERROR

8. HOW TO PERFORM SPIROMETRY EXAMINATION

Put the spirometry attachment cable into the patient cable socket. Switch on the unit with the use of the  key. The display shows:

STOP	I	II	III
25mm/s	10mm/mV		

To obtain spirometry option press the  key. The display shows:

	SEX	AGE	HEIGHT
PATIENT:	>M	40	168

Set the patient data with the use of "←", "→" keys. Set ">" mark in required position and set:

sex (*Male* or *Female*), age, height with the use of „↓“ „↑“ keys.

Press the  key. The display shows:

EXPIRE	SEX	AGE	HEIGHT
PATIENT:	>M	40	168

Make expiration - the display shows:

WAIT	SEX	AGE	HEIGHT
PATIENT:	>M	40	168

The printout is made.

It is possible to carry out a copy of the spiro examination if it is any stored in the unit memory.

To obtain the copy of spirometry examination press the  key in the spiro option. To return to ECG mode simply press the  key.

9. THE UNIT WITH ALPHANUMERIC KEYBOARD - MODEL B, D.

9.1 How to insert the patient data.

The unit equipped with alphanumeric keyboard has the possibility to insert the patient data such as: name, surname, age, height, sex.

1. To insert the patient data:
in STOP mode press **Enter**. The display shows:

SURNAME

It is possible to Insert the following patient data:

- surname - name
- ID - age
- height - weight
- sex - blood pressure
- physician
- 1* (place for comments:
- 2* two lines)

To change the parameters press **Enter**. To leave from the mode press **Esc**.

2. To obtain the numbers press **Shift** and **required number** at the same time.
3. To accept press **Enter** "↓".

4. To exit press **Esc** "↑" or any functional key.

Use:

Enter to insert the mode

Esc to exit the mode

"←" to delete the last letter

Ctrl and "←" to delete one parameter (for example "Surname")

Ctrl and **Esc** to delete all of the patient data parameters

10. HOW TO MAINTAIN THE UNIT.

To clean the cover use a duster moistened with water and soap.

DO NOT USE SOLVENT AND CORROSIVE LIQUID ! ! !

11. ELECTROCARDIOGRAPH AsCARD A4 EQUIPMENT

1. Limb electrodes (EKK)* 4 pcs
2. Chest electrodes (EPP)* 6 pcs
3. ECG Patient cable - 1pc
4. Main cable - 1 pc
5. Paper 112 mm - 1 roll
6. ECG gel - 250 ml
7. Operation manual - 1 pc

* the electrodes (type EKK and EPP) and other accessories are always available at Aspel.

12. TYPICAL PROBLEMS OF THE UNIT.

The unit doesn't switch on	- there is no power supply - the internal accumulator is load out Switch the unit on to power supply - check the switcher on/off - check the safety device
The unit stops during the registration	- the internal accumulator is load out Switch the unit on to power supply – the internal accumulator will load. It is possible to make registration on internal accumulator
The unit signalises INOP	Any of electrode can be disconnected. Check the electrodes. If there is still INOP the patient cable can be damaged. Please contact Authorised Service.
ERROR STROB	The unit was kept in too low temperature. Wait until the temperature increases.
The printout in AUTO mode impossible	There was incorrect data inserted to the unit memory. Delete RAM memory – switch off the unit with the use of  key, then switch on the unit with use of  and  key pressed at the same time.
SPIRO – the printout independent on the patient's exhaust	There was incorrect patient data inserted. Delete RAM memory – switch off the unit with the use of  key, then switch on the unit with use of  and  key pressed at the same time
OUT OF PAPER	Put the new roll of paper to the paper box of the unit.
	The accumulator is discharged. Charge the accumulator

TABLE OF CONTENTS

1..	APPLICATION	4
1.1	The description of the unit.	5
1.2	Picture of the unit	5
1.3	Basic Technical Parameters.	6
1.4	Manufacturer and Representative.	7
1.5	Version of the unit.	8
2.	WORKING CONDITIONS.	9
3.	FEATURES.	9
3.1	How to prepare AsCARD A4 to work.	10
3.2	How to put ECG paper roll to the unit.	11
3.3	The arrangement of the keys.	12
3.4	The arrangement of electrodes.	14
4.	EXAMINATION PERFORMANCE.	16
4.1	How to perform manual examination.	16
4.2	How to perform automatic examination.	17
4.3	How to perform copy of automatic examination.	19
5.	HOW TO CHANGE REGISTRATION PARAMETERS.	20
5.1	Record speed.	20
5.2	Record sensitivity.	20
5.3	ECG leads.	20
5.4	Additional filters.	21
5.5	Registration mode.	21
6.	HOW TO CONFIGURATE THE UNIT PARAMETERS.	22
7.	HOW TO TRANSMIT ECG EXAMINATION BY FAX.	26
8.	HOW TO PERFORM SPIROMETRY EXAMINATION.	27
9.	THE UNIT WITH AN ALPHANUMERIC KEYBOARD – MODEL B	28
9.1	How to insert the patient data.	28
10.	HOW TO MAINTAIN THE UNIT.	29
11.	ELECTROCARDIOGRAPH AsCARD A4 EQUIPMENT.	30
12.	TYPICAL PROBLEMS WITH THE UNIT.	31